

RELAY[®] PRO


THORACIC STENT-GRAFT SYSTEM

Uniquely Inspired
for Ideal Placement



Achieving Precision with the Relay®Pro Delivery System for Accurate and Controlled Deployment

Relay®Pro employs a delivery system engineered to achieve a **perpendicular** positioning of the stent graft at the proximal landing zone, designed for **optimal apposition** at the inner aortic curve.



“Relay®Pro’s ability to land accurately combined with its low profile will allow me to successfully treat complex anatomy with precision.”¹



WATCH
Relay®Pro
Deployment



WATCH
Relay®Pro NBS
Deployment

- 1 Proximal Clasping**
 - ▶ **Repositioning** of the stent-graft for precise placement
 - ▶ **Facilitating** accurate and perpendicular deployment
- 2 Precurved inner catheter**
Conforms to the aortic arch designed for alignment of the stent-graft
- 3 Inner Sheath**
Allows for releasing the stent’s energy in phases, for more accurate placement of the stent-graft
- 4 Controller**
Allows for staged deployment enhancing control and accuracy in stent-graft placement
- 5 Mechanical Advantage**
Forward and backward gear system allows for small incremental movements of the stent-graft enhancing controlled delivery

Empowering Confidence in Deploying NBS Configuration


Relay®Pro NBS delivery system implements **two features** that have been designed for precise and safe proximal deployment, **minimizing birdbeak** and **retroflex** effects

Two support wires* guide the inferior portion toward the inner aortic wall, keeping it aligned with the landing zone, **minimizing the risk of retroflex**

* NOTE: The support wires are only present in stent-grafts with proximal diameters of 32mm or greater.



The **Flared End** configuration of the inner sheath is designed for proper alignment and to **minimize birdbeaking**



“Accurate deployment with favorable apposition even in hostile aortic arches contributed to low rates of early and mid-term complications.”^{2 ^}



Preoperative CTA with 3D reconstruction illustrates the benefit of a low-profile delivery system⁴

100%
Technical
success^{*1,3}
110/110
56/56

Thoracic Aortic Aneurysm and Penetrating Atherosclerotic Ulcer Cohort
Acute Complicated Type B Aortic Dissection Cohort

0%
Birdbeak through
12 months^{**3}
0/56

Acute Complicated Type B Aortic Dissection Cohort

** These studies include all Relay®Pro with the NBS configuration being predominant



Postoperative CTA 3D reconstruction shows RelayPro NBS conformability to an angulated aortic arch⁴

[^] The Relay®Pro indications for use will vary by region. Note that in the US and Canada, Relay®Pro is only indicated for treatment of the DTA and not specifically indicated for IMH. Always consult IFU.

² El Beyrouiti *et al.* (2020). Early results of a low-profile stent-graft for thoracic endovascular aortic repair. *PLOS ONE*

³ Rossi *et al.* (2024). One-Year Results of a Low-Profile Endograft in Acute, Complicated Type B Aortic Dissection. *The Annals of Thoracic Surgery*

⁴ Case images courtesy of Prof. Wilson Szeto, Cardiovascular Surgery, Penn Presbyterian Medical Center, <https://www.vumedi.com/video/2nd-stage-tevaring-with-thoraflextm-hybrid-staying-on-label-at-all-times/>



Discover solutions for every segment of the aorta
[terumoaortic.com](https://www.terumoaortic.com)



LinkedIn



VuMedi



X