

RELAY[®] PRO


THORACIC STENT-GRAFT SYSTEM

Uniquely Inspired
for Ideal Placement



Achieving Precision with the Relay®Pro Delivery System for Accurate and Controlled Deployment

Relay®Pro employs a delivery system engineered to achieve a **perpendicular** positioning of the stent graft at the proximal landing zone, designed for **optimal apposition** at the inner aortic curve.



“Relay®Pro’s ability to land accurately combined with its low profile will allow me to successfully treat complex anatomy with precision.”¹



WATCH
Relay®Pro
Deployment



WATCH
Relay®Pro NBS
Deployment

1 Proximal Clasp

- **Repositioning** of the stent-graft for precise placement
- **Facilitating** accurate and perpendicular deployment

2 Precurved inner catheter

Conforms to the aortic arch designed for alignment of the stent-graft

3 Inner Sheath

Allows for releasing the stent’s energy in phases, for more accurate placement of the stent-graft

4 Controller

Allows for staged deployment enhancing control and accuracy in stent-graft placement

5 Mechanical Advantage

Forward and backward gear system allows for small incremental movements of the stent-graft enhancing controlled delivery

Empowering Confidence in Deploying NBS Configuration


Relay®Pro NBS delivery system implements **two features** that have been designed for precise and safe proximal deployment, **minimising birdbeak** and **retroflex** effects

Two support wires* guide the inferior portion toward the inner aortic wall, keeping it aligned with the landing zone, **minimising the risk of retroflex**

* NOTE: The support wires are only present in stent-grafts with proximal diameters of 32mm or greater.



The **Flared End** configuration of the inner sheath is designed for proper alignment and to **minimise birdbeaking**



100%

Accurate device deployment^{** ^ 2}

23/23

Aortic dissection, aortic aneurysms, PAUs and IMH; N=1 aortic erosion and aortic rupture, each. The RelayPro is NOT indicated for erosion or rupture.

0%

Birdbeak through 12 months^{** 3}

0/56


Acute, Complicated Type B Aortic Dissection Cohort.

** These studies include all Relay®Pro with the NBS configuration being predominant


“Accurate deployment with favorable apposition even in hostile aortic arches contributed to low rates of early and mid-term complications.”^{2 ^}




86 yo, PAU/focal dissection in the DTA⁴



Final Angio after Relay®Pro NBS implant⁴



Post Operative CT-SCAN⁴



[^] The Relay®Pro indications for use will vary by region. Note that in the US and Canada, Relay®Pro is only indicated for treatment of the DTA and not specifically indicated for IMH. Always consult IFU.

² El Beyrouti et al. (2020). Early results of a low-profile stent-graft for thoracic endovascular aortic repair. *PLOS ONE*

³ Rossi et al. (2024). One-Year Results of a Low-Profile Endograft in Acute, Complicated Type B Aortic Dissection. *The Annals of Thoracic Surgery*

⁴ Case images courtesy of Wilson Y. Szeto, Chief, Division of Cardiovascular Surgery, Hospital of the University of Pennsylvania-Penn Presbyterian, <https://www.vumedi.com/video/relaypro-thoracic-stent-graft-features-to-clinical-practice/>



Discover solutions for every segment of the aorta
[terumo-aortic.com](https://www.terumo-aortic.com)



LinkedIn



VuMedi



X